

5304530069 FREEZER THERMOMETER, EAST WEST MANUFACTURING (SKU# L304432837)

5304530069 THERMOMÈTRE DE CONGÉLATEUR, EAST WEST MANUFACTURING (UGS L304432837)

5304530069 TERMÓMETRO PARA CONGELADOR, EAST WEST MANUFACTURING (SKU N.º L304432837)

Installation

1. Place the thermometer as close to the center of the freezer as possible.
2. Do not allow the thermometer to touch anything but the freezer shelf on which it is placed.
3. Allow 24 hours before reading the temperature.
4. The temperature will rise every time the door is opened.

Installation

1. Placez le thermomètre aussi près que possible du centre du congélateur.
2. Ne laissez pas le thermomètre toucher quoi que ce soit d'autre que la tablette du congélateur sur laquelle il est placé.
3. Attendez 24 heures avant de prendre une lecture de la température.
4. La température augmente chaque fois que la porte est ouverte

Instalación

1. Coloque el termómetro lo más cerca posible del centro del congelador.
2. Deje que el termómetro toque solo el estante del congelador en el que se coloca.
3. Espere 24 horas antes de leer la temperatura.
4. La temperatura aumentará cada vez que se abra la puerta.